



BS-23-11-27-001-V2.2

Scaling Check




L ►

M ►

S ►

Full Face Mask

F5+ / F5A+ / F5 / F5A / F5AS /
F4 / HM01 / F3 / BMC-FM2/
F2 NV1 / F2 NV2 / F1B / BMC-FM1A



- 01 Align the edge with the dent of chin
- 02 Check aligned position with nasal bridge to determine the size

family cares




BS-23-11-27-001-V2.2

Scaling Check



Nasal Mask

N5+ / N5A+ / N5B+ / N5 / N5H /
N5A / N5AH / N5B / N5BH / NM4





BS-23-11-27-001-V2.2

Scaling Check



Nasal Pillows Mask

P6 / P2 / P2H

