



Exhaust-flow information

F&P Evora™ Full Mask No diffuser		F&P Vitera™ Mask No diffuser		F&P Simplus™ Mask No diffuser		F&P Evora Nasal Mask No diffuser		F&P Eson™ 2 Mask Diffuser		F&P Eson 2 Mask No diffuser		F&P Solo™ Nasal Mask No diffuser		F&P Nova Micro™ Nasal Pillows Mask No diffuser		F&P Solo Nasal Pillows Mask No diffuser		F&P Brevida™ Mask Diffuser		F&P Brevida Mask No diffuser	
Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)	Pressure (cmH ₂ O)	Flow (L/min)
4	21	4	24	4	22	4	17	4	19	4	27	4	15	4	13	4	15	4	21	4	25
6	25	5	26	5	25	5	19	5	22	5	31	6	20	6	16	6	20	5	24	5	27
8	28	6	28	6	27	6	21	6	25	6	34	8	25	8	19	8	25	6	27	6	30
10	32	7	30	7	29	7	23	7	27	7	36	10	29	10	22	10	29	7	29	7	32
12	35	8	32	8	31	8	24	8	29	8	39	12	34	12	24	12	34	8	31	8	34
14	38	9	34	9	33	9	26	9	31	9	41	14	38	14	27	14	38	9	33	9	37
16	41	10	36	10	35	10	28	10	34	10	44	16	42	16	29	16	42	10	35	10	39
18	44	11	38	11	37	11	29	11	36	11	46	18	46	18	31	18	46	11	37	11	41
20	46	12	40	12	39	12	31	12	38	12	48	20	49	20	33	20	49	12	39	12	42
22	49	13	42	13	40	13	32	13	39	13	50	22	49	22	35	22	49	13	41	13	44
24	51	14	43	14	42	14	34	14	41	14	52	24	49	24	37	24	49	14	43	14	46
26	54	15	45	15	43	15	35	15	43	15	54	26	49	26	38	26	49	15	44	15	48
28	56	16	47	16	45	16	37	16	45	16	56	28	49	28	40	28	49	16	46	16	49
30	58	17	48	17	46	17	38	17	47	17	57	30	49	30	42	30	49	17	47	17	51
		18	50	18	48	18	39	18	48	18	59			32	44	32	49	18	49	18	52
		19	51	19	49	19	40	19	50	19	61			34	46	34	49	19	50	19	54
		20	53	20	50	20	41	20	51	20	62			36	48	36	50	20	52	20	55

F&P Brevida operates up to 25 cm H₂O.
To view the complete exhaust flow table refer to the F&P Brevida Use & Care Guide.
F&P Eson 2 operates up to 30 cm H₂O.
To view the complete exhaust flow table refer to the F&P Eson 2 Use & Care Guide.

F&P Evora operates up to 25 cm H₂O.
To view the complete exhaust flow table refer to the F&P Evora Use & Care Guide.
F&P Nova Micro operates up to 25 cm H₂O.
To view the complete exhaust flow table refer to the F&P Nova Micro Instructions for Use.

F&P Simplus operates up to 25 cm H₂O.
To view the complete exhaust flow table refer to the F&P Simplus Use & Care Guide.
F&P Solo operates up to 20 cm H₂O.
To view the complete exhaust flow table refer to the F&P Solo Use & Care Guide.

F&P Vitera operates up to 30 cm H₂O.
To view the complete exhaust flow table refer to the F&P Vitera Use & Care Guide.

F&P, Brevida, Eson, Evora, Nova Micro, F&P Solo, Simplus and Vitera are trademarks of Fisher & Paykel Healthcare Limited. For patent information see www.fphcare.com/ip.
SUI 621512 REV D 2023-11 © 2023 Fisher & Paykel Healthcare Limited